


Wireless Probe Type Ultrasound Scanner

Welcome to inspect our probe with your phone/ipad

GG
ایده گستر قدرت



wifi wireless



A new generation of magic medical devices

Features:

- Workable with Tablet PC or Smart Phone (iOS, Android, Windows)
- Built-in and replaceable battery
- Advanced digital imaging technology, clear image
- High cost-effective
- Wireless connectivity, easy to operate
- Small and light, easy to carry
- Applicable in hospital, emergency, clinic, outdoor and vet inspection
- Intelligent terminal platform, powerful expansion functions on application, storage

Standard Configuration:

- Wireless probe device 1 Unit
- Built-in Battery and charge cable 1 Set
- User Manual copy

Specifications:

- Scan mode: Electronic Array, Electronic Convex array scanning;
- Frequency: convex probe 3.5MHz
- Elements: 128 elements
- Scanning Depth: Convex 100-305mm, adjustable
- Field of view: Convex R60
- Screen: Smart phone or tablet screen
- Supporting system: iOS, Android, Windows
- Display mode: B, B/M, Color, PDI, PW
- Frame rate: 18f/s
- Image gray scale: 256 level
- Image/Video Storage: Store on mobile phones, Tablet PC
- Measure: distance, area, obstetrics and other;
- Power: by built-in lithium battery, 4200mAh
- Battery working time: 3.5 hours
- Size: 156mm×60mm×20mm
- Net Weight: 250 gram



Sonostar Technologies Co., Limited



Address: C Building, #27 Yayingshi Road, Science Town, Guangzhou, China www.sonostar.net
Tel:+86-20-32382095 Fax:+86-20-62614030 Skype:sonostar-ut Mail:sales@sonostar.net