


Wireless Probe Type Ultrasound Scanner

Welcome to inspect our probe with your phone/ipad



A new generation of magic medical devices



Scan to get
detail introduction



Features:

- Small and small size, easy to carry.
- Wireless type without probe cable, working freely.
- Waterproof design, convenient for sterilization.
- Remote diagnostic convenience, capability and easy image transfer

Standard Configuration:

- Wireless probe device: 1 Unit
- Built-in Battery and charge cable 1 Set
- User Manual copy

Specifications:

- Scan mode: Micro-convex or linear array scanning;
- Frequency: Micro-Convex probe 5.0/7.0MHz, Linear probe 10/12MHz
- Elements: 128 elements
- Scanning Depth: Micro-Convex 60-160mm, Linear 20-80mm, adjustable
- Field of view: Micro-Convex R20, Linear L25
- Screen: Smart phone or tablet screen
- Supporting system: iOS, Android, Windows
- Display mode: B, B/M, Color, PDI, PW

- Frame rate: 18f/s
- Gain: 30db-105db
- Dynamic Range: 40/50/60/70/80/90/100/110
- Noise reduction: 0-1-2-3-4
- Image gray scale: 256 level
- Image/Video Storage: Store on mobile phones, Tablet PC
- Measure: distance, area, obstetrics and other;
- Power: by built-in lithium battery, 2800mAh
- Battery working time: 2.5hours
- Size/Net Weight: 156mm×60mm×20mm, 250g



Sonostar Technologies Co., Limited

Address: C Building, #27 Yayingshi Road, Science Town, Guangzhou, China
Tel: +86-20-32382095 Fax: +86-20-62614030 Skype: sonostar-ut



www.sonostarmed.com
Mail: sales@sonostar.net