

IGG

ایده گستر قدرت



More information, please browse www.sonostarvet.com

Sonostar Technologies Co., Limited

Factory: First floor, Investment Building, Nanping, Zhuhai, Guangdong, China
Market Center: 504#, C Building, #27 Yayingshi Road, Science Town, Guangzhou
Tel: +86-20-32382095 Fax: +86-20-62614030 Skype: sonostar-ut
www.sonostarvet.com Mail: sales@sonostar.net MSN: sonostar@163.com

Note: Sonostar trademarks is protected by Madrid Treaty.

SonoStar™


Professional Create Value



Light, clear image!

V5

Palm Ultrasound Scanner

- 7 inch TFT screen
- Pretty and very light
- Full digital imaging technology, crystal-clear Image
- Can work with backfat probe
- Low power consumption
- 2G Large volume storage image
- Built-in battery
- Solid shell and structure
- Waterproof(Optional)



Waterproof(optional)



Varius Probes



Backfat Probe

Standard Configuration:

Host, convex probe, built-in battery, power adapter, bag.

Optional:

Linear probe, Micro-convex probe, Rectal probe, Backfat probe.



Specification:

- Display screen: 7" TFT
- Scanning manner: Electronic convex array, electronic linear array
- Working frequency: 2.0MHz~10MHz
- Display mode: B, B+B, B+M, M, 4B
- Scanning range: Convex array 60°~150°
- Amplification factor 1.0, 1.2, 1.5, 2.0
- Gain control: 8 segments TGC and overall gain can be adjusted respectively
- Image processing: 8 γ corrections, frame correlation, point correlation, line correlation, digital filtering, digital edge enhancement and pseudo color processing, etc
- Conventional measurements: distance, perimeter, area, volume
- Obstetric measurement: horses, cattle, sheep, dogs, cats obstetric measurement
- Image gray-scale: 256 levels
- Scanning line number: 512 lines/frame
- Frame rate: 30 frames/second
- Digital scanning conversion: 512×512×8 bits
- Battery capacity: 3000mAh/7.4V, can work more than 2 hours
- Memory capacity: ≥2GB
- Report function: auto make report
- Output interface: USB port, video.
- Overall size: 216mm×147mm×27mm
- Weight: 950g(with battery)
- Waterproof: IP65(optional)



Full digital imaging technology bring crystal-clear image