

Wireless Probe Type Ultrasound Scanner

Welcome to inspect our probe with your phone/ipad



L5TC



IGG
یده گستر قدرت



Specifications:

- Probe Type: Consist by two of Convex, Transvaginal
- Scanning mode: Electronic array
- Display mode: B, B/M, color doppler version with B+Color, B+PDI, B+PW
- Probe element: 128
- Channel of RF circuit board: 32
- Probe frequency and scan depth, head radius/width, scan angle(convex):
Convex head 3.5MHz/5MHz, 90/160/220/305mm, 60mm, 60°
Transvaginal head 6.5MHz/8MHz, 40/60/80/100mm, 13mm, 149°
- Image Adjust: BGain, TGC, DYN, Focus, Depth, Harmonic, Denoise, Color Gain, Steer, PRF
- Cineplay: auto and manual, frames can set as 100/200/500/1000
- Puncture assist function: the function of in-plane puncture guide line, out-of-plane puncture guide line, automatic blood vessel measurement.
- Measure: Length, Area, Angle, heart rate, Obstetrics
- Image save: jpg, avi and DICOM format
- Image frame rate: 18 frames/second
- Battery working time: 2.5 hours(according to whether keep scan)
- Battery charge: by USB charge or wireless charge, take 2 hours
- Dimension: L156×W60×H270mm
- Weight: 250g
- Wifi type: 802.11g/20MHz/5G/450Mbps
- Working system: Apple iOS and Android, Windows



Sonostar Technologies Co., Limited FDA TUV CE 0123
 Address: C Building, #27 Yayingshi Road, Science Town, Guangzhou, China www.sonostar.net
 Tel:+86-20-32382095 Fax:+86-20-62614030 Skype:sonostar-ut Mail:sales@sonostar.net

WWW.IGG-MED.COM

شرکت ایده گستر قدرت نمایندده انحصاری فروش و خدمات پس از فروش سونوگرافی موبایل کمپنی سونوگستر ایران