



Wireless Probe Type Ultrasound Scanner

Welcome to inspect our probe with your phone/ipad

شرکت ایده گستر قدرت توانمند انحصاری فروش و خدمات پس از فروش گیمینی SONOSTAR



L5PC



A new generation of magic medical devices

Features:

- Advanced digital imaging technology, clear image
- High cost-effective
- Small and small size, easy to carry.
- Wireless type without probe cable, working freely.
- Waterproof design, convenient for sterilization.
- Remote diagnostic convenience, capability and easy image transfer

Standard Configuration:

- Wireless probe device 1 Unit
- Built-in Battery and charge cable 1 Set
- User Manual copy

Specifications:

- Scan mode: Electronic Array, Electronic Linear array scanning;
- Frequency: Linear probe 10/14MHz
- Scanning Depth: Linear 20/30/40/55mm, adjustable
- Screen: Smart phone or tablet screen
- Display mode: B, B/M, Color, PDI, PW
- Image gray scale: 256 level
- Image/Video Storage: Store on mobile phones, Tablet PC
- Measure: distance, area, obstetrics and other;
- Power: by built-in lithium battery, 4200mAh
- Size: 156mm×60mm×20mm
- Elements: 128 elements
- Field of view: Linear L25
- Supporting system: iOS, Android, Windows
- Frame rate: 18f/s
- Battery working time: 3 hours
- Net Weight: 220 gram


ایده گستر قدرت



تراسنکین



Needle puncture guide



Android tablet/iPad



Tablet bracket



Trolley



Sonostar Technologies Co., Limited



Address: C Building, #27 Yayingshi Road, Science Town, Guangzhou, China www.sonostar.net
Tel:+86-20-32382095 Fax:+86-20-62614030 Skype:sonostar-ut Mail:sales@sonostar.net